

PTFE Loop Prosthesis for Stapedotomy and Stapedectomy

CODE
SPL 03.16/17/18/19

PTFE: Polytetrafluorethylene

Conventional PTFE prosthesis featuring major technical developments of the hook.

Seat for guided opening:
open the hook slightly and press the seat into the long process of the incus.
The seat prevents a too wide opening of the hook which could compromise the elastic memory of the PTFE.

Seat for forceps:
the seat shows the exact point of hold and allows for a firm grip.

Sized Prostheses

Total length (L)	Piston diameter			
	Ø 0.5 mm	Ø 0.6 mm	Ø 0.4 mm	Ø 0.8 mm
4.00 mm	SPL 03.16.400	SPL 03.17.400	SPL 03.18.400	SPL 03.19.400
4.50 mm	SPL 03.16.450	SPL 03.17.450	SPL 03.18.450	SPL 03.19.450
5.00 mm	SPL 03.16.500	SPL 03.17.500	SPL 03.18.500	SPL 03.19.500
5.50 mm	SPL 03.16.550	SPL 03.17.550	SPL 03.18.550	SPL 03.19.550
6.00 mm	SPL 03.16.600	SPL 03.17.600	SPL 03.18.600	SPL 03.19.600
6.50 mm	SPL 03.16.650	SPL 03.17.650	SPL 03.18.650	SPL 03.19.650
7.00 mm	SPL 03.16.700	SPL 03.17.700	SPL 03.18.700	SPL 03.19.700


Trimvable Prostheses

Dimensions

Total length: 7.0 mm
Trimvable piston diameter:
0.4 - 0.5 - 0.6 - 0.8 mm

Packaging

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Piston diameter mm	Code
	7.00	7.00	0.4
0.5			SPL 03.16
0.6			SPL 03.17
0.8			SPL 03.19

PTFE Loop Prosthesis for Stapedotomy and Stapedectomy

CODE
SPL03.16x/17x/18x/19x

PTFE: Polytetrafluorethylene

Conventional PTFE prosthesis featuring major technical developments of the hook.

Seat for guided opening:
open the hook slightly and press the seat into the long process of the incus.
The seat prevents a too wide opening of the hook
which could compromise the elastic memory of the PTFE.

Sized Prostheses

Total length (L)	Piston diameter			
	Ø 0.5 mm	Ø 0.6 mm	Ø 0.4 mm	Ø 0.8 mm
4.00 mm	SPL 03.16.400x	SPL 03.17.400x	SPL 03.18.400x	SPL 03.19.400x
4.50 mm	SPL 03.16.450x	SPL 03.17.450x	SPL 03.18.450x	SPL 03.19.450x
5.00 mm	SPL 03.16.500x	SPL 03.17.500x	SPL 03.18.500x	SPL 03.19.500x
5.50 mm	SPL 03.16.550x	SPL 03.17.550x	SPL 03.18.550x	SPL 03.19.550x
6.00 mm	SPL 03.16.600x	SPL 03.17.600x	SPL 03.18.600x	SPL 03.19.600x
6.50 mm	SPL 03.16.650x	SPL 03.17.650x	SPL 03.18.650x	SPL 03.19.650x
7.00 mm	SPL 03.16.700x	SPL 03.17.700x	SPL 03.18.700x	SPL 03.19.700x


Trimvable Prostheses

Dimensions

Total length: 7.0 mm
Trimvable piston diameter:
0.4 - 0.5 - 0.6 - 0.8 mm

Packaging

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Piston diameter mm	Code
	7.00	0.4	SPL 03.18x
0.5	SPL 03.16x		
0.6	SPL 03.17x		
0.8	SPL 03.19x		